

GREENCARE A1C

CLINICAL CHEMISTRY: WET CHEMISTRY CATEGORY

BEST VALUE FOR MONEY

One Step HbA1c measuring system

PROVIDE ACCURATE QUANTITATIVE TEST RESULTS FOR HBA1C
NGSP and IFCC certified and CV is less than 3%.

PROVIDE THE MOST CONVENIENT SOLUTION.
Store and operate the kit at room temperature.

Help clinician decision during the patient's visit.



TOUCH SCREEN



BUILT-IN
THERMAL PRINTER



BUILT-IN
THERMAL PRINTER



COMMUNICATE
WITH LIS



VOICE GUIDE

MUREX Diagnostic Products Specialists



ISO 9001:2015
CERTIFIED

08 Bonifacio Street, Poblacion East,
San Nicolas 2447, Pangasinan &
119 14th Ave., Cubao, Quezon City, Metro

+639 189647158
+075-600-8159
+02-8-962-8743

@murexdiagnostics.ph
www.murexdiagnostics.com
info@murexdiagnostics.com



SCAN ME



SCAN ME



Provide accurate and precise results / Certified NGSP

Accuracy

Model. RT 100

Present high accuracy and excellent correlation compared with HPLC with 40 different samples provided by NGSP.

GRENCARE A1c System was evaluated with 40 different samples provided by NGSP. The results show high accuracy and excellent correlation.



Difference, %	GRENCARE A1c System
within±3%	50%(20/40)
within±6%	92.5%(37/40)
within±10%	100%(40/40)

Precision

Model. RT 100


Within-Run Precision (IFCC sample)
Each of the samples was measured 20 times per each for precision using 2 meters across 3 lots.

IFCC calibrator	Normal	Abnormal
Mean HbA1c%	5.5	8.9
CV (%)	2.4	1.9

Day-to-day Precision (IFCC sample)
Each of the samples was measured 1 time a day for 20 days using 2 meters across 3 lots.

IFCC calibrator	Normal	Abnormal
Mean HbA1c%	5.4	9.0
CV (%)	2.6	1.9


Simple Operation without preparation



1 Insert the cassette.



2 Collect the blood from the reagent pack.



3 Insert the reagent pack and the result is released in less than 5 minutes.

GreenCare A1c Specification

GreenCare A1cAnalyzer	
Model	RT-100
Dimension	205(W) × 200(L) × 150(H)mm
Weight	1,600g
Power	100/240VAC 12VDC, 3.33A
Display	4.3" LCD Resolution(480 X 272 pixel)
Memory	400 tests
Printer	Internal thermal printer
Additional function	USB x2, Ethernet, LIS connected External barcode reader available
Operating condition	Temperature: 15 ~ 35°C / Humidity: 15 ~ 75%
Storage condition	Temperature: -10 ~ 60°C / Humidity: 15 ~ 75%

GreenCare A1c Test Kit	
Measurement item	Hemoglobin A1c
Principal	Boronate affinity chromatography
Accuracy	Within ± 6% 92.5%
Precision	< 3% CV
Measurement Range	4.0 ~ 15.0% / 20 ~ 140 mmol/mol
Sample type	Capillary whole blood Venous blood with anticoagulant (EDTA, Heparin)
Sample Volume	3.5µl
Measuring time	Less than 5 minutes
Storage condition	2 ~ 30
Expiration date	12 months from manufacture

MUREX Diagnostic Products Specialists



ISO 9001:2015
CERTIFIED

08 Bonifacio Street, Poblacion East,
San Nicolas 2447, Pangasinan &
119 14th Ave., Cubao, Quezon City, Metro

+639 189647158
+075-600-8159
+02-8-962-8743

@murexdiagnostics.ph
www.murexdiagnostics.com
info@murexdiagnostics.com



SCAN ME



SCAN ME