

URIT *Since 1984*
Global Diagnostics Supplier



UC-58

Semi-automatic Urine Analyzer

- Up to 550 T/H
- 7-inch extra 50% resolution sensitive touch screen
- Available from 1 to 14 parameters including CR, MA, CA, VitC
- Real strip image display
- Temperature control system integrated
- Lithium battery enhanced



URIT Medical Electronic Co., Ltd.

Head Office
Add.: No. D-07 Information Industry
District, High-tech Zone, Guilin, Guangxi, 541004, P.R. China
Tel: +86-773-2288586 Fax: +86-773-2288560
E-mail: export@uritest.com service@uritest.com
<http://www.urit.com>

Shenzhen Office
Add.: Room EJ, 24th Floor, New Baohui Building, Nanshan Road,
Nanshan District, Shenzhen, P.R. China
Tel: +86-755-26413721

Rev: 202101-1 (20210118)



Guilin, China



UC-58

Semi-automatic Urine Analyzer



Flexible Test Mode

- Continuous measuring 550 T/H
- Strip automatic unloading



14 Items Strips

- 14 test items including CR, MA, CA, VitC
- Albumin-to-Creatinine Ratio (ACR) helps to screen early kidney disease
- Color and turbidity results integrated



Complete QC System

- URIT QC solution provides all test items, including VitC



Advanced CIS System

- High accuracy and high efficiency
- 5 wavelengths guarantee precise and confident results
- Real image display



Powerful

- DC 12V & Lithium battery integrated
- More than 600 tests available without power supply



Easy to Use

- Touchscreen with linux system
- Analyzer automatically detects strip existence, type and starts tests



Constant Temperature Control System

- Avoid temperature interference
- Allow accurate testing in minimum to 5°C room temperature
- Offer the best reaction environment

Parameters 18 Parameters Leukocyte, Ketone, Nitrite, Urobilinogen, Bilirubin, Protein, Glucose, Specific Gravity, Blood, pH, Ascorbic Acid, Microalbumin, Creatinine, Calcium, ACR, PCR, Color, Turbidity

Measurement Principle Reflectance photometry

Reaction Time 45 seconds

Throughput Up to 550 T/H

Display 7-inch high resolution (1024*600) capacitive color touch screen

Memory Capacity 20,000 results

Strip HC series strips

Interface Serial port (Rs232), USB (barcode reader), wifi, bluetooth optional

Work Condition Environment temperature: 5~30°C, RH≤80%, lumination≤2000lx

Power Supply Instrument: 12V DC; Adapter: 100V-240V, 50/60Hz

Dimension & Weight 308 mm (L) × 210 mm (W) × 151 mm (H), 3 kg